



Technical resources

Agate grinder

Model: MNY

<https://www.tencanballmills.com/en/products/grinding-series/other-grinding-equipment/agate-grinder.html>

09.10.2026 / TENCAN

Agate grinder

Model: MNY

Ideal lab equipment for grinding small samples. Agate mortar ensures contamination-free fine powder. Widely used in research labs.

Technical specifications

Mortar opening diameter	Φ120mm
Mortar speed	10rpm
Pestle speed	20rpm-90rpm

These specifications describe the series. Confirm the selected model, power supply and configuration before ordering.



The agate mortar grinder is mainly composed of three parts: a grinding rod, an agate mortar and a scraper. During the working process, the agate mortar rotates clockwise and the grinding rod crushes the solid particles in the mortar to make the solid particles smaller and finer while being fixed. The scraper scrapes off the particles on the edge of the mortar, so that the solid particles return to the mortar and continue grinding. During the movement of the agate mortar and grinding rod, the grinding opportunity for each solid particle is equal, so the particle size of the ground micron powder is uniform. Set a reasonable grinding time for different solid particles (generally, the longer the grinding time, the finer the particles are ground), so that both coarse particles and fine particles are ground into very fine powder particles. The particle size can reach the micron level, and some can even reach the nanometer level. The entire grinding process can be automatically controlled. Agate bowl mill is an ideal grinding equipment for laboratory grinding of small amounts of micron powder.

研磨系列

产品细节



Technical parameters

model	MNY-Φ120
Voltage	AC110V/220V 50Hz/60Hz
Grinding rod power	10w
Mortar power	25w
Mortar and pestle material	agate
Mortar diameter	Φ120mm
Mortar depth	34mm
Mortar speed	10rpm
Rotation speed of grinding rod	20rpm-90rpm
Grinding time	0-900min

1. The power of grinding particles can reach micron level.
2. Precisely control the grinding rod speed and grinding time.
3. Used for grinding superhard material powder, electronic powder, resin powder, food powder, pharmaceutical powder, and chemical powder.
4. Simple operation.

Contact

Changsha Tianchuang Powder Technology Co., Ltd. / Changsha, Hunan, China

Email: lichang0516@gmail.com

Phone / WhatsApp: +86 731 84027560 / WhatsApp +86 19021571516

<https://www.tencanballmills.com/en/products/grinding-series/other-grinding-equipment/agate-grinder.html>

Before operation, review the selected model's instructions and the material requirements.